



**DUFNER®**  
Instruments



Robotic Surgery  
Bulldog Clamps

## Features

- ▶ Easy access through single port (SP), flexible trocar (FT) or rigid trocar (RT)
- ▶ Worldwide unique Bulldog clamps which can be placed, manipulated and removed by the robot and the endo applicator
- ▶ Two modes: distal or lateral placement
- ▶ Closing pressure: ~170 grams
- ▶ Single surgeon placement with Pro-Grasp® Instrument



**20919-35-353**

DeBakey Bulldog clamp, 50 mm serration, curved, 102 mm total length, access through SP or FT without sleeve



**20919-35-354**

DeBakey Bulldog clamp, 58 mm serration, straight, 110 mm total length, access through SP, FT, RT



**20919-35-355**

DeBakey Bulldog clamp, 30 mm serration, straight, 83 mm total length, access through SP, FT, RT



**20919-35-356**

DeBakey Bulldog clamp, 30 mm serration, curved, 83 mm total length, access through SP, FT, RT



**20919-35-363**

DeBakey Bulldog clamp, 56 mm serration, slightly curved, 109 mm total length, access through SP, FT, RT



**20919-35-364**

DeBakey Bulldog clamp, 50 mm serration, Satinsky, 102 mm total length, access through SP or FT



**20919-35-365**

DeBakey Bulldog clamp, 70 mm serration, Satinsky, 125 mm total length, access through SP or FT



**26932-60**

MICASEPT Applicator for robotic Bulldog clamp, for lateral pickup, pistol grip with ratchet, Ø 10 mm, 350 mm working length



**26932-62**

MICASEPT Applicator for robotic Bulldog clamp, for distal pickup, pistol grip with ratchet, Ø 10 mm, 380 mm working length




**26910-13**

Trocar set, pyramidal tip with folding valve, with insufflation tap, Ø 13 mm, 100 mm working length



**20913-17**

MICTEC Flexible trocar set TITANIUM, Ø 13 mm, 70 mm long



DUFNER Instrumente GmbH  
Föhrenstraße 9  
D-78532 Tuttlingen  
Phone: +49 (0) 7461/3697  
Mail: [info@dufner-tuttlingen.de](mailto:info@dufner-tuttlingen.de)  
[www.dufner-tuttlingen.de](http://www.dufner-tuttlingen.de)